

HEMATOLOGY CONTROLS



- ▶ Ready to use, stable and homogenous controls made from fresh human blood like patients' samples
- ▶ Ensure quality patient results
- ▶ Reliable & assess analyzer accuracy
- ▶ Precision in quantitative blood analysis
- ▶ Available in Low, Normal and High levels

Volume	2ml/3ml	Shelf Life	105 days
Open Vial Stability	21 days	Storage Temperature	2-8° C

AVAILABLE FOR



DYNACOUNT 3D

Automated **3 Part**
Differential Hematology
Analyzer

DYNACOUNT 5D

Automated **5 Part**
Differential Hematology
Analyzer



CPC Diagnostics Pvt Ltd

4th floor, Westminster, 108, Dr. Radhakrishnan Salai,
Mylapore, Chennai - 600 004

Phone : +91 9840462339, +91 8754468400

Email : support@cpcdiagnostics.in